

## *Determination of preservatives in cosmetic products*

**Test aim:** to determine the presence of the 12 preservatives in matrices of cosmetic origin (i.e.: make up, body and tooth care products, sun products, skin creams and emulsions, etc...).

**Procedure:** the extraction of the preservative agent is performed in acid environment through an ethanol/water mixture. The substance is heated at 60°C for 5 minutes and cooled (4-5 °C) in fridge for an hour in order to facilitate the precipitation of possibly present fatty acids. After filtration, preservatives will be determined through HPLC with UV-VIS detector at variable wave length (Diode Array). The following table shows the various preservatives according to their elution and corresponding detection limit and range.

<b>Preservative</b>	<b>Detection limit (mg/kg)</b>	<b>Range (mg/kg)</b>
<i>4-Hydroxybenzoic acid</i>	10	10-5000
<i>Salicylic acid</i>	50	50-20000
<i>Benzoic acid</i>	100	100-20000
<i>2-phenoxyethanol</i>	100	100-20000
<i>Sorbic acid</i>	10	10-5000
<i>Methyl Paraben</i>	10	10-5000
<i>Dehydroacetic Acid</i>	20	20-5000
<i>1-phenoxypropan-2-ol</i>	50	50-20000
<i>Ethyl-Paraben</i>	10	10-5000
<i>Propyl-Paraben</i>	10	10-5000
<i>Butyl-Paraben</i>	10	10-5000
<i>Benzyl-Paraben</i>	10	10-2500

**Instruments:** HPLC - DAD

**Method:** MP-0806 R2 internal standard